

Sleeve print proof v6

TRUE SIZE • BUILT TO VITALL'S ARTWORK SPEC • 23 AUGUST 2026

The v5 concepts transposed onto the real dimensions taken from Vitall's own artwork file (**2025 Box Design.eps**, InDesign, authored 14 Aug 2025). Pages 2 to 4 are the three kit sleeves as flats, at true size, in the orientation a printer receives them. Page 5 is the back face shown upright so the compliance window can be checked against a ruler.

DIMENSIONS

Flat trim	179.00 × 304.00 mm — from the file's own HiResBoundingBox (507.4016 × 861.7323 pt)
Folded face	179.00 × 152.00 mm , two faces, lower one rotated 180°
Fold	Horizontal centre line, 152.00 mm from trim top. Assumed , see below
Bleed	3 mm all round → page is 185 × 310 mm
Safe type area	5 mm inside trim

THE COMPLIANCE WINDOW – LIFTED FROM THE VECTOR, NOT MEASURED

Size	108.00 × 57.00 mm
From left trim edge	34.09 mm
From right trim edge	36.91 mm — not centred , sits 2.8 mm left of centre
Vertical	Centred on its face, 47.50 mm clear above and below
Share of panel	22.7%
Source	5-point closed path at x 96.63–402.77 pt, y 135.35–296.93 pt

WHAT CHANGED FROM V5, AND WHY

- **The window grew by 27.0 mm in width and moved 19.5 mm down.** v5 reserved 81.0 × 51.5 mm at 49.0 mm from the left; the real rectangle is 108.00 × 57.00 mm at 34.09 mm. That is 47% more area.
- **v5's activation block collided with the real window.** It sat at 11–32 mm across and 90–111 mm down, and the window occupies 47.5–104.5 mm down. The back face is re-laid-out into three bands: brand above the window, window clear, regulatory below.
- **The QR code is gone from the sleeve back.** Per the /activate deprecation already recorded in 02_brand/STATE.md: the sampling QR belongs on the **insert**, and "activate your kit" becomes "see your results", because the kit is already linked to the customer at order creation.
- **Fronts are unchanged** apart from being scaled to true size. Only the sleeve front is kit-specific; the back is universal across all three.

TWO THINGS THIS PROOF ASSUMES

1. The fold is at the midpoint. Vitall's file carries no cut or crease layer — no spot colours, no custom colours, one page — so it is printed artwork, not a cutter guide. 152.00 mm is exactly half of 304.00 mm and the geometry supports it, but it has not been confirmed.

2. Only the back face carries a window. A rectangle of identical size sits on the other face in Vitall's artwork. If **both** faces need one, the hero front has a 108 × 57 mm hole in the middle of it and pages 2 to 4 need redrawing, not adjusting. **This is the first question to put to Ben.**

ANDROPRIME



AT-HOME HORMONE TEST

Testosterone

ANALYSED IN A UKAS-ACCREDITED UK LABORATORY

Thresholds signed off by Dr Ewa Lindo, GMC-registered GP

IVD CE • UKCA REF • AP-T01

ANDRO-PRIME.COM

FOLD (assumed) • 152.00

ANDRO-PRIME.COM

In-vitro diagnostic for self-collection. Read the enclosed instruction sheet before use.
Distributed by Andro Prime Ltd (Co. 17185839), UK. • IVD • CE • UKCA
Analysed in a UKAS-accredited UK laboratory.
Store below 25°C. Do not use if the pouch is damaged.



36.91

34.09

SEE YOUR RESULTS AT ANDRO-PRIME.COM
Your kit is already linked to you. Nothing to register.
Read the enclosed instruction sheet before taking your sample.

ANDROPRIME

ANDROPRIME



AT-HOME ENERGY & RECOVERY TEST

Energy & Recovery

ANALYSED IN A UKAS-ACCREDITED UK LABORATORY

Thresholds signed off by Dr Ewa Lindo, GMC-registered GP

IVD CE · UKCA REF · AP-E02
ANDRO-PRIME.COM

FOLD (assumed) · 152.00

ANDRO-PRIME.COM

In-vitro diagnostic for self-collection. Read the enclosed instruction sheet before use.
Distributed by Andro Prime Ltd (Co. 17185839), UK. · IVD · CE · UKCA
Analysed in a UKAS-accredited UK laboratory.
Store below 25°C. Do not use if the pouch is damaged.



SEE YOUR RESULTS AT ANDRO-PRIME.COM
Your kit is already linked to you. Nothing to register.
Read the enclosed instruction sheet before taking your sample.

ANDROPRIME

ANDROPRIME



AT-HOME HORMONE & RECOVERY TEST

Hormone & Recovery

ANALYSED IN A UKAS-ACCREDITED UK LABORATORY

Thresholds signed off by Dr Ewa Lindo, GMC-registered GP

IVD

CE · UKCA

REF · AP-H03

ANDRO-PRIME.COM

FOLD (assumed) · 152.00

ANDRO-PRIME.COM

In-vitro diagnostic for self-collection. Read the enclosed instruction sheet before use.
Distributed by Andro Prime Ltd (Co. 17185839), UK. · IVD · CE · UKCA
Analysed in a UKAS-accredited UK laboratory.
Store below 25°C. Do not use if the pouch is damaged.



Read the enclosed instruction sheet before taking your sample.


Your kit is already linked to you. Nothing to register.

SEE YOUR RESULTS AT ANDRO-PRIME.COM

ANDROPRIME

Sleeve back — reading orientation

UNIVERSAL ACROSS ALL THREE KITS • TRUE SIZE • MEASURE THIS PAGE TO VERIFY



Band layout, forced by the window

- Above 0 – 47.50 mm brand + what to do next
- Window 47.50 – 104.50 mm left clear for Vitall's lot / expiry / barcode / kit ID
- Below 104.50 – 152.00 mm regulatory fine print

Two deliberate changes from v5

- The QR is gone. Per the /activate deprecation the sampling QR belongs on the insert.
- "Activate your kit" is now "see your results" — the kit is linked at order creation.